



Spot Installation Instructions for Windows

This document explains the procedure for manual installation of Spot for Windows platforms. The installation procedure is as follows:

1. If R is not installed, fetch it from the following website and install it.
<http://cran.r-project.org/bin/windows/>
2. If imview is not installed, fetch it from the following website and install it.
http://www.hca-vision.com/Spot_Documentation/imview-1.0.1-win32.exe
3. Copy the Spot package to your hard disk. You may have downloaded it from the following website.
http://www.hca-vision.com/product_download.html
4. Start R and select the "Packages->Install Package from local zip file" option. Use the browse option to select the Spot_3.1.zip file.
5. Exit R, log out of Windows and log back in again.
6. Start R and execute the following commands from inside R – “library(Spot)” and “runSpots()”. The Spot package should be loaded and display a license message, and a graphic user interface should appear.